





Inclusion Support Program

Specialist Equipment Catalogue

Contact us:

- t 1800 129 606
- e inclusion@gowriesa.org.au

Gowrie SA Inclusion Agency 42 Crittenden Rd, Findon SA 5023





INDEX

Seating	5
Tables, Trays	11
Standing Frames	12
Toileting Equipment	14
Hoists and Slings	17
Bathing Support	18
Ramps and Accessories, Swings	19
Switch / Cause and Effect Toys, etc.	20
A Box Full of Feelings	22
How to Access Specialist Equipment	23







Terms and Conditions

DISCLAIMER: The information in this catalogue is intended to provide general information on Specialist Equipment available to loan to eligible Early Chilhood Education and Care (ECEC) services. This catalogue can be used as a starting point for Therapists and ECEC educators in gaining knowledge about resources available from Gowrie SA.

Please note:

- At date of printing, information in this Catalogue is accurate, however product specifications and availability may change without prior notice.
- This is not an extensive Equipment list. Please contact us with enquiries about additional items.
- Please also note Equipment received may not be the ones shown in the images, however they will operate in the same manner as the requested Equipment.

Specialist Equipment is an equipment that facilitates and supports the inclusion of a child (or children) with additional needs in the care environment. Gowrie SA considers an equipment being 'Specialist', when it is valued over \$200. The eligible ECEC services are responsible for the purchase of any equipment that is commercially available and valued below this amount.

All eligible ECEC services are able to borrow Specialist Equipment, where a need is identified by a Therapist or Inclusion Professional (IP).

Borrowing most Specialist Equipment requires a written recommendation by a Therapist (e.g. Occupational Therapist or Physiotherapist, etc.) A limited number of non-adjustable and less advanced Specialist Equipment can be recommended by an IP, without a therapist's written recommendation.

Specialist Equipment may be borrowed by emailing inclusion@gowriesa.org.au or phone 1800 129 606.

All eligible ECEC services accessing Specialist Equipment must submit a Request for Service and have a Strategic Inclusion Plan (SIP) in which the need for Specialist Equipment must be identified.

Borrowers must fill in a Specialist Equipment Request Form and sign a Loan Agreement Form which covers all loans. The loan period will depend on the needs of the child, the service, and the type of equipment provided. For finding the latest version of this catalogue and our request forms please use this link: http://www.gowriesa.org.au/request-forms

All items borrowed must be used under adult supervision. Please ensure you adhere to safety guidelines when using Equipment.

For all posture-related equipment and many of more advanced items, a Therapist must adjust the Equipment for the child at least once in 12 months, or as the need arises.

Gowrie SA, as the Inclusion Agency (IA) for South Australia will arrange the delivery of Specialist Equipment within South Australia.

Specialist Equipment can take some time when it is on order so please contact us as soon as you are made aware of a need for it.

For equipment repairs and modification, please contact Gowrie SA as soon as possible.



Seating



Kelly Chairs

Ideal for home, schools and preschools use with children who need some support when sitting.

The chair can be adapted with a large range of accessories. (Abduction Wedge, Footplate, 'H' Harness, Pelvic Strap, Mobile Base, Adjustable Headrest, Lateral Pads, Leg Skids, Anti Tilt Stabiliser Bars, etc.)

- Pre School Kelly Chair (K2) (Pictured)
- Primary School Kelly Chair (K3)



CAP Chairs

These chairs can be customised with adjustable footplate, armrests, pommel and pelvic and shoulder straps.

Other accessories include: Trays, Dual Locking Rubber Castors and Plastic Castors, Sliders, Thoracic supports, Head Rest extensions.

Pre-school and Primary sizes available.



Jenx Corner Chair

Has an adjustable height back and can be used with or without the upper back support. An abduction block is used to help keep the user well positioned, with their hips placed well back into the seat.

Also includes adjustable table (photo under trays/tables, page 10)

Size 1 – 6 months to 2 years old, Max weight: 25kg

Size 2 – 1 year to 6 years old, Max weight: 35kg

Size 3 – 4 years to 10 years old, Max weight: 50kg





Tumble Form Floor Sitters with Base and Tray

Includes the Feeder Seat Positioner and the Floor Sitter Wedge.

The Wedge can be used to keep the Positioner in an upright or a reclined position suitable for the child's

A Velcro brand hook & loop system holds the seat securely, allowing instant adjustability.

Available with adjustable Tray in Small (User Height: up to 91cm) and Medium (User Height: up to 122cm).





Special Tomato Soft Touch Sitter

Floor sitters that are Latex Free, seam free, peel-resistant and impermeable to fluids.

Design provides up to 25° of tilt. Adjustable 5-point harness. Contoured head & lateral support.

Size 1: Target Age: 1-2 years; Target Height: 76-89cm Size 2:Target Age: 2-6 years; Target Height: 89-122cm Size 3: Target Age: 6-9 years; Target Height: 101-142cm



The Tumble Forms Ready Racer is recommended for children 18 months to five years. The Ready Racer is a pre-wheelchair, paediatric mobility aid that combines high-performance wheelchair handling with sports-car styling.

Inner seat dimensions: 24 x 58 x 34 cm Weight limit: 25 kg







Leckey PAL Classroom Seats (left)

Complexity: Mild

Suitability and Dimensions:

Size 1: 1-3 years, max 30 kg; Seat Depth: 18-24cm;

Seat Width: 15-25cm

Size 2: 2-5 years, max 35 kg; Seat Depth: 22-29cm;

Seat Width: 17-27cm

Size 3: 4-8 years, max 40 kg; Seat Depth: 25-33cm;

Seat Width: 23-33cm





Smirthwaite Juni Chair

For children who require minimal levels of postural support when seated. Comes with contoured seat and backrest, built-in lumbar support, pelvic strap, height and depth adjustable seat, height adjustable armrests as well as angle adjustable backrest as standard. (Back Angle: 95/100/105). Also look at the Smirthwaite Juni Table (page 11).

Size 1: H: 16 - 22 cm; W: 24 cm; D: 18 -24 cm; Back Height: 26 cm Size 2: H: 20-28 cm; W: 28 cm; D: 22-30 cm; Back Height: 29 cm Size 3: H: 24-34 cm; W: 32.5 cm; D: 26-36 cm; Back Height: 32.5 cm

Kidoo Chair

Positioning chair with adjustable seat to backrest angle, fixed tilt in space adjustment, adjustable seat depth & chair height.). Also look at the Kidoo Table (page 11).

Size 1: W: 23 cm; D: 18-23.5 cm; Back Height: 26 cm Size 2: W: 26.5 cm; D: 22.5-30 cm; Back Height: 31 cm Size 3: W: 31 cm; D: 24.5-34.5 cm; Back Height: 35 cm

Size 4: W: 35 cm; D: 31-41 cm; Back Height: 40cm







Jenx ATOM Chair on Height Adjustable Y Base

Max user weight: 25 kg Age range: 6 Months to 4 Years Seat depth: 9-34cm

Seat width: 18-27cm

Seat height from floor: 33.5-75.5cm Lower leg length: 16.5-32.5cm

Back Pad height: 27-38cm

Chest width: 20-28cm

Back angle: 15° - 15° (posterior - anterior)

Tilt in space angle: 30° - 10° (posterior - anterior)





Jenx Bee Chair on Height Adjustable Y Base

Bee is a seating system specifically for younger children.

Age Range (Approx.): 4 months - 5 Years Back Angle (Supine-Prone):10 Degrees Back Pad Height (Seat to Top): 27 - 36 cm

Chest Width: 15 - 27 cm

Depth with Long Leg Conversion: 37 - 56 cm

Hip Width: 13 - 25 cm Lower Leg Length: 9.5 - 34.5 cm

Seat Depth: 11 - 35 cm

Tilt-in-Space Angle (Supine-Prone): 25 Degrees

Max user weight: 25 kg

Y-base Size (Footprint): 65 - 78 cm

Seat Height from Floor (Y-base): 28.5 - 75 cm



Jenx Gamma Chair on Y Base

Size 1: 1-5 years; Size 2: 3-8 years

Max weight 45kg

Specifications:

Seat to floor: 33.5-52 cm

Seat length 18.5-26.5cm (size 1); 22-32cm (size 2)

Tilt-in-space (prone/recline): 15°Prone - 25°Recline

Back column angle (prone/recline): -15° to +15°

Armrest to floor: 41.5-76.5cm (size 1); 44-79cm (size 2)


Thoracic pad width: 17-26cm

Hip width: 16-33cm (size 1); 18-33cm (size 2)

Lower leg length: 17-31cm (size 1); 17-36cm (size 2)

Back height: 35-45.5cm (size 1); 39-49.5 (size 2)









Rifton Activity Chair, Small hi/lo base (R830)

Tool free adjustments. Wide range of optional accessories available to create an adaptive chair providing customised active seating:

Frame width: 66 cm

Seat height above floor: 25-64 cm

Seat angle (tilt-in-space): 15° forward, 25° back

Backrest angle: 5° forward, 20° back Seat height above foot-board: 23–30 cm Seat width with hip guides: 18-23 cm

Seat depth: 20-30 cm

Armrest height above seat: 13–19 cm Trunk support width: 14–29 cm Backrest height: 12½–15½ cm

Headrest height above seat: 37–53 cm

Max. working load: 34 kg



Jenx Junior Seating system on Height Adjustable Y Base

Jenx Junior is a seating system for mild, moderate & complex needs. Size 1 (18 months—5 years); Size 2 (3–8 years)

Features:

Outstanding Pelvic Control System;

Leg Control Pads (controlled independently of the pelvic position);

PU (polyurethane leather) Support Pads;

Adjustable Seat Depth;

Back Angle Adjustment (15 degrees recline through to 15 degrees prone);

Adjustable Back Height (31-41 cm for size 1; 33-43 cm for size 2);

Max User Weight: 50kg



Leckey Mygo Seating System Size 1

Modular Mobility Seating

Complexity: Moderate to Complex

Age (approx.): 3 - 10 years User Weight: 18-50kg User Height: 105-150cm Seat Width: 20-32.5cm Seat Depth: 27-42cm Knee Width: 9-11cm

Footplate: Abduction: 8°; Adduction: 12°;

Plantarflexion/Dorsiflexion: 10°

Chest Width: 17-27cm Backrest Height: 36-47cm

Backrest Angle: Prone: 10°; Recline: 25°

Seat to Sandal: 21.5-35cm Top of Seat to Floor: 36-70cm







Special Tomato Height Right size 1

3-Point Pelvic Belt Included Easily adjustable seat plate and foot plate. Made from natural hardwood Wooden Tray with protective cover Soft-Touch Seat and Back Cushions added for additional comfort





Zoomi Plus High Chair

Height range: 11 – 61.5cm

Width: 45cm Depth: 55cm Back height: 26cm

Back angle: 95°/100°/105°

Frame height: 81

Maximum user weight: 85kg



Leckey Squiggles Early Activity System

Age: 0-36 Months

Complexity: Mild - Complex

Designed to be used as part of a programme of early intervention

to maximise a child's future potential.

Consists of a number of soft, bright rolls, wedges, straps and supports, which when used in combination with the Activity Development Programme by therapists and parents, enhance the development of babies and young infants.



Kid'in'Form Positioning System

Kit includes:

comfort mat, cylindrical cushions (2 small, 2 medium), lateral positioning cushions, belt straps (Medium and Large), backrest, large cylindrical roll, pierced cushion



Tables and Trays



CAP Corner Seat Table Small (above left) and Large (above right)

Has rails to suit corner chair table.

Measurements for Small: (Could be used for Jenx Corner Chairs Size 1 & 2)

L: 45cm; W: 70cm; H: 16 - 25cm

Measurements for Large: (Can be used with Jenx Corner Chair Size 3) Depth: 50cm; Width 80cm; Height Range 30-50cm, with 5cm increments



Jenx Corner Table

This table is fully adjustable and helps to encourage a good posture. It also offers an ideal surface for a range of play and educational activities. Can be used with the Jenx Corner chair or stand alone.

Table Top Angle Range (From Horizontal): 40° - 40°

Size 0: Table Size: 51x60.5cm; Height 19-28cm; pitch between legs: 57cm Size 1: Work area: 49x30cm; Height 24.5-34.5cm; pitch between legs: 57cm Size 2: Work area: 59x28cm; Height 24.5-40cm; pitch between legs: 67cm Size 3: Width: 70cm; Depth: 50cm; Height Range: 30-45cm; Cut-out dimen-

sions: D: 16cm: W: 32cm

Smirthwaite Juni Table Size 1

The legs are adjustable so the child can work or play at the correct height and the table folds, which allows for easy storage.

Height: 31 - 45 cm Width/Length: 70 cm

Depth: 53 cm





Kidoo Table Size 1

W: 65cm; D: 43cm; H: 39 - 50cm

Tilt angle: 0° – 12°

Maximum user weight: 30kg





Standing Frames





Dynamic Standing Frame

Adjustable heights, pelvic and leg support Maximum Working Load: 60kg Working Height: 45-60cm (Small); 63-83cm (Medium) Vertical adjustable knee strap; Non-slip Rubber Mat; Pelvic Support: Variable supportive pads providing total pelvic support

Base Dimensions: 43 x 59cm

Leckey Totstander Size 1

Adjustable height, pelvic and leg support Complexity: Mild to Moderate Suitability: Age 1-5 years Technical Overview: Base: Depth and Width: 45cm Maximum User Weight: 22kg





Jenx Penguin Upright Stander

A flexible and adjustable standing system, providing assisted upright positioning for younger children with moderate ability.

Build in table tray with bowl for sensory activities Age Range (approx.): 18 months - 4 years

Base Size: 60 x 50cm

Chest Height (footplate to top of chest pad): 60-80cm

Chest & Hip Pad Width: 14-25cm

Max User Weight: 25kg

Tray Dimensions: Depth: 36cm; Width: 53cm

Jenx Kangaroo Upright Stander

A flexible and adjustable standing system, providing assisted upright positioning for older children with moderate ability. Has the capacity to adjust to offer true upright standing. Has adjustment in three planes: knees, pelvis and chest.

Build in table tray with bowl for sensory activities Age Range (approx.): 3-7 year

Base Size: 70 x 58.5cm

Chest Height (footplate to top of chest pad): 70-100cm

Chest Pad Width: 20.5-30cm Hip Pad Width: 21-30cm

Max User Weight: 35kg

Tray Dimensions: Depth: 39cm; Width: 58cm





Jenx Monkey Prone Standing Frame

A highly adjustable standing system that provides Prone to

Upright positioning

Build in table tray with bowl for sensory activities

Angle Range: 20°- 90°

Age Range (approx.): 9 months - 4 years

Base Size: 57.5 x 88cm

Chest Height (footplate to top of chest pad):

51-75cm

Chest & Hip Pad Width: 15-23.5cm

Max User Weight: 25kg

Tray Dimensions: Depth: 36cm; Width: 53cm





Jenx Pronestander size 1

Offers a continuous range of angles so the optimum position for each child can be achieved. The individual support pads can be adjusted in relation to each to ensure every child can be properly supported.

Build in table tray with bowl for sensory activities

Age Range (approx.): 3-7 years

Base Size: 90x 65cm

Tilt-in-Space (supine) (90°=upright): 40°-0°

Max User Weight: 55kg

Chest Width & Hip Width: 19-29cm

Hip Pad Height: 41-75cm

Jenx Multistander: Supine/Prone/Upright Standing System Size 1

Provides highly supportive standing in prone, supine or upright. Beneficial to growing children in multi-use environments. Has a wide variety of head, trunk, leg and foot support options. The deliberately open design ensures you can always see how the user is positioned.

Built in table tray

Angle Range: 10° -90°

Age Range (approx.): 9 months-6 years

Max User Height: 116cm

Frame Dimensions: 56.5cm (W) x 71cm (L)

Chest Pad Width: 15-26cm

Hip Pad Width: 12-26cm

Height of Stander (floor to top of upper pad): 60-90cm







Toileting Equipment



Toilet Frame with Step

Adjustable Rail height, Available in different step heights 14-20cm. Tube style frame with step. Fits over most toilets.





Adjustable Height Toilet Frame Surround

Rails:

Minimum height: 55cm (approximately)
Maximum height: 70cm (approximately)

Width: 55cm

Depth: 70cm from front of unit to rear of unit



Lite Toilet Step - Adjustable Legs (CAP)

With or without castors (right)
Standard step size: 48 cm wide x 60 cm deep
Minimum height: 10.5 cm with 4.5 cm of height adjustment
Minimum height, using the longer legs: approx. 15 to 20 cm
Rail height: 48 cm to 80 cm from top of platform
Distance between rails: 41 cm
Handrails adjust in 3 cm increments from 50-80 cm high.
Lockable Rubber Wheeled Castors on Toilet Side of Step,
Fits most standard pedestals





Smirthwaite Combi adjustable toileting chair Size 1 and 2

With standard seat, potty and toilet seat with bucket, hip belt, adjustable seat depth.

Could be used as a supportive chair whilst providing numerous features to assist with the everyday toileting needs. Seat height and depth can be adjusted to suit the needs of child.





Nuvo Petite Children's Toilet Platform

Overall Height: 54.7 – 66.7cm

Overall Width: 56cm Overall Depth: 57.2cm

Platform Height: 9.6 – 21.6cm Aperture Width: 22.5cm Aperture Depth: 21.5cm User Weight Limit: 76kg



A lightweight, portable toilet seat. For ages 2-9 years.

One frame for both seat sizes; Facilitates a supported squat posture, making toileting easier.









Electric Change Tables

Adjustable height for comfortable and ergonomic use Electrical change table with lockable castors and adjustable side rails. Padded vinyl top.

- · ABCO Electric Change Table- Junior (1.4m & 1.2m with side cots)
- Alrick Electric Change Table (1.2m, 1.4m & 1.8m)
- · ArjoHuntleigh Electric Change Table- Omega (1.8m)

Folding Portable Athlegen Change Table

Height adjustable change table. Two sizes available.

W: 70 cm L: 135 cm or 180 cm

H: 60-75cm or 65-80cm









Padded and Plastic Toilet Seat Reducer Ring

For a slender or young child who may require extended toileting time. Offers comfort and non-slip security. Has a large abductor/splash-guard that is padded and removable. Fits any toilet size. This Model has a small hole size: 18 x 20cm.



Throne rail - Child Standard Size

Toilet rail + 2 bolts with rubber caps Handle height: 19 cm Between handles: 40 cm Height above seat: 28 cm

Child back plate internal measurements: 29 cm

Also available: Throne Grab Rail 3 in 1



Toilet Steps With Non-Slip Top

Plywood toilet step: varnished; non slip top surface $56 \times 56 \text{ cm}$ 10 or 16 cm high, with 20 x 10 cm cut-out to fit toilet



Hoists and Slings



Hoist Lifters and Accessories

Mobile hoist, rechargeable battery operated. Safe working load capacity from 147kg - 155kg.

- · Hoist Oxford Advance Patient Lifter Folding
- · Hoist Battery Charger
- · Hoist Lifter Battery (Removable)
- · Hoist/Lifter Wangara W147
- · Kerry Equipment Home Care Lifter (KH402)

Assorted Slings for Hoists

Tech Assist - Mesh Full Body With Head Support- Small, X Small & XXS

· Commode Independent - Large

Regent Full Body X-Small, Small & Medium

Europa Full Body Mesh Sling For Hoist- Small

All Purpose Mesh - Extra Extra Small

Junior Full Body Mesh Europa- Large

Invacare General Purpose- With Head Support - Small & XS

· General Purpose Small- Mesh

Mesh Full Body With Head Support







Bathing Support



Flamingo Shower Chair / Toilet Seat Size 2

A toilet/bathing chair. Can be tilted and the back can be angled and therefore assist children with special needs with their seating posture.

The Flamingo is easy to clean and adjustable in height. It can be placed over the toilet or can be used with a pan.

Angling: Seat: -15° to 35° Back: -90° to 45° Foot


support bar: 90° to 60°

Seat width: 28cm Seat depth: 28.5cm Back height: 37.5cm Lower leg length: 23-33.5 Seat height: 49.5-59.5cm Back width: 21.5cm Width: 45.5cm Length, frame: 66cm Height frame w/seat: 89cm Max. user weight/load: 50kg





Ramps and Accessories



DECPAC Ramps

(Various sizes, widths and lengths)

- · Decpac 90 cm
- · Decpac Folding 135 cm Edge Barrier Limiter
- · Folding Ramp with Gutters 165 cm









Switch/Cause and Effect Toys





Switches - in different sizes and colours . Limited switch toys available. (below and next Page) Includes animals: Cassy Bear; Percy the Pig; Robbie the Rabbit; Peter Rabbit; Dave the Dalmatian; Pip the Penguin; Lucy the Dog, etc.





Bright Red Switch



Lighted Musical Tunes



Carnival Tube Switch







Kaleido Dome with Shapes





Bead Chain in a Box



Big Mack





Talkable 4





A box full of feelings

Designed specifically to support the social-emotional development of children aged 2-7.
Built around four basic feelings: happiness, fear, anger and sadness.

Contains several basic materials to get children involved in more than 20 different activities.





How to access Specialist Equipment

https://gowriesa.org.au/inclusion-agency/specialist-equipment-library

- · Specialist Equipment is free to borrow for eligible services, including delivery.
- · Borrowing most Specialist Equipment requires a written recommendation by a Therapist (e.g. Occupational Therapist or Physiotherapist, etc.) A limited number of non-adjustable and less advanced Specialist Equipment can be recommended by an IP, without a therapist's written recommendation.
- · All eligible ECEC services accessing Specialist Equipment must submit a Request for Service and have a Strategic Inclusion Plan (SIP) in which the need for Specialist Equipment must be identified.
- Borrowers must fill in a Specialist Equipment Request Form and sign a Loan Agreement Form which covers all loans. The loan period will depend on the needs of the child, the service, and the type of equipment provided. For finding the latest version of this catalogue and our request forms please use this link: http://www.gowriesa.org.au/request-forms
- For all posture-related equipment and many of the more advanced items, a Therapist must adjust the Equipment for the child before use and at least once in 12 months, or as the need arises.
- For equipment repairs and modification, please contact Gowrie SA as soon as possible.







Gowrie SA **Inclusion Agency**

- t 088234 5219 (Option 4) e inclusion@gowriesa.org.au
- 42 Crittenden Road, Findon SA 5023

